



BeneSphera™ Brand

ELISA PLATE WASHER

W21



Features and Benefits

- Flat, U-bottom or V-bottom plates and strip washing
- Liquid level warning
- Emergency stop available
- Includes an 8-way manifold; 12-way manifold available as an option
- Low residual volume
- Large LCD display with Win CE compatible operating system
- Large memory to store up to 50 user-programmed wash protocols



System Specifications

Resolution	50µl
Dispense precision	<5% CV at 350µl; +-2% at 300µl
Residual Volume	<3µl (standard); <1µl(double aspiration)
Display	90 mm X 53 mm LCD
Data Storage Capacity	50 user-programmed wash protocols
Manifold	8 pins standard
Operating Condition	10°C - 40°C
Store Condition	1°C - 40°C
Dimensions	405X343X175 mm
Weight	12 Kg
Power Supply	AC 200-240V, 50/60Hz/150W max



Phillipsburg, NJ 9001:2008 & 14001:2004
Paris, KY 9001:2008
Mexico City, Mexico 9001:2008
Deventer, the Netherlands 9001:2008, 14001:2004 & 13485:2003

Gliwice, Poland 9001:2008 & 17025:2005
Selangor, Malaysia 9001:2008
Dehradun, India 9001:2008, 14001:2004 & 13485:2003
Mumbai, India 9001:2008 & 17025:2005



About Avantor™ Performance Materials

Avantor Performance Materials manufactures and markets high-performance chemistries and materials around the world under several respected brand names, including the J.T.Baker®, Macron Fine Chemicals™, Rankem™, BeneSphera™ and POCH™ brands.



Avantor products are used in a wide range of industries. Our biomedical and life science solutions are used in academic, industry and quality control laboratories for research, pharmaceutical production and medical lab testing, while our electronics solutions are used in the manufacturing of semiconductors and flat panel displays. Based in Center Valley, Pennsylvania (USA), Avantor is owned by an affiliate of New Mountain Capital, LLC.

For additional information please visit www.avantormaterials.com or follow [www.twitter.com/avantor_news](https://twitter.com/avantor_news)

Ordering Information and Assistance

Customer Service and Technical Service

India

1201-1206, 12th Floor, Pinnacle Business Tower,
Shooting Range Road, Surajkund, Faridabad-121009, India

TEL: +91-129-4267000, 4267100, 4267200

FAX: +91-129-4267299, 4267199

E-MAIL: diagnostics.india@avantormaterials.com

SERVICE TOLL FREE NO: 1800-102-5555

Netherlands

Teugseweg 20, 7418 AM Deventer, The Netherlands

TEL: +31-570-687500

FAX: +31-570-687574

E-MAIL: avantor.emea@avantormaterials.com

www.avantormaterials.com

ASK Avantor™

Our Web site features ASK Avantor™, which includes live chat capabilities with customer service representatives.
www.avantormaterials.com/askavantor

Corporate Headquarters

Avantor Performance Materials, Inc.
3477 Corporate Parkway
Suite #200
Center Valley, PA 18034
USA

Worldwide Locations

• China • Malaysia • North America
• India • Mexico • Poland
• Korea • The Netherlands • Taiwan

For contact information at these locations, visit
www.avantormaterials.com/WorldwideDirectory

Lit # 9129

©2013 Avantor Performance Materials, Inc. All rights reserved.

Trademarks are owned by Avantor Performance Materials, Inc. or its affiliates unless otherwise noted.