



ID-Centrifuge 12 S II



ID-Centrifuge 12 S II

Product Specifications

- Centrifuge head for 12 ID-Cards
- Monitored by microprocessor
- Automatic balance control of the centrifuge head
- Rpm and time are displayed (LCD) in English, French, German, Italian, Spanish or Portuguese
- Pre-fixed centrifugation speed and time at 1,030 rpm \pm 5 rpm (85 g \pm 1 g) for 10 minutes
- Dimensions: 30 x 36 x 18 cm (W x D x H)
- Weight: 6.2 kg
- Electrical data: 110-240 V/50-60 Hz
- CE compliant, according to IVD Directive 98/79/EC

Ordering Information

Catalog No.	Description	
009550	ID-Centrifuge 12 S II (complete with head)1 unit
009551	Accessory head1 unit

The Complete Solution for Safe Transfusion





Bio-Rad Laboratories is a multinational manufacturer and distributor of clinical diagnostics and life science research products.

*Clinical
Diagnostics Group*

DiaMed
1785 Cressier FR
Switzerland
www.bio-rad.com/immunohematology • Click on ID-System