



ID-Incubator 37 S I



ID-Incubator 37 S I

Product Specifications

- For the incubation of ID-Cards at 37°C
- Capacity of 2 ID-Card racks (24 ID-Cards)
- Monitored by microprocessor
- Standard program: 15 minutes at 37°C (± 2°C)
- Programmable incubation time from 1-60 minutes
- Time and temperature are displayed (LCD) in English, French, German, Italian, Spanish or Portuguese
- Dimensions: 31 x 34 x 19 cm (W x D x H)
- Weight: 6 kg
- Electrical data: 110-240 V/50-60 Hz
- CE compliant, according to IVD Directive 98/79/EC

Ordering Information

Catalog No.	Description
009690	ID-Incubator 37 S I1 unit

The Complete Solution for Safe Transfusion





Bio-Rad Laboratories is a multinational manufacturer and distributor of clinical diagnostics and life science research products.

*Clinical
Diagnostics Group*

DiaMed
1785 Cressier FR
Switzerland
www.bio-rad.com/immunohematology • Click on ID-System