

Thanks to your input Bio-Rad has developed its next generation of centrifuges. ID-Centrifuge L is the latest addition to the range and has been designed to meet your specific needs. The main advantage is that you can use both ID-Cards and tubes in a single instrument.

Our centrifuge technology helps you to standardize your procedures. All our centrifuges are CE compliant according to IVD Directive 98/79/EC.

ID-Centrifuge L

This instrument is the latest centrifuge developed by Bio-Rad which can be used with both ID-Cards and tubes.

Technical Data

Performance


RCF maximum	1,000 g (for tubes)
RCF	85 g (for ID-Cards)

Dimensions

Width	380 mm
Height	230 mm
Depth	500 mm
Weight	13 kg
(heads not included)	

Main Features

- Only one centrifuge for tubes and ID-Cards
- Automatic detection of the head type
- Automatic selection of the appropriate program
- Head for post transportation ID-Card centrifugation with adapted program (5 minutes at 85 g)
- 28 tube capacity
- Panel PC touch-screen with easy-to-use interface (icon based)
- 4 compatible centrifuge heads
- Multi-language software
- CE compliant according to IVD Directive 98/79/EC



Ordering Information

Catalog #	Description
ID-Centrifuge L	
009201	ID-Centrifuge L, 1 unit
009410	Head 24 ID-Cards, 1 unit
009412	Head 28 tubes, 1 unit
009413	Head 12 tubes, 1 unit
009414	Head 2 racks, 1 unit

Thanks to your input Bio-Rad has developed its next generation of incubators. ID-Incubator L is the latest addition to the range and has been designed to meet your specific needs. The main advantage is that you can use both ID-Cards and tubes in a single instrument.

Our incubator technology provides you with security in the incubation process and standardizes your manual work. All our incubators are CE compliant according to IVD Directive 98/79/EC.

ID-Incubator L

This instrument is the latest incubator developed by Bio-Rad which can be used with both ID-Cards and tubes.

Technical Data

Electrical data

Voltage 110-240 VAC
Frequency 50/60 Hz

Dimensions

Width 380 mm
Height 230 mm
Depth 500 mm
Weight 13 kg
(heads not included)

Main Features

- Incubator for both tubes and ID-Cards
- Independent incubation areas (ID-Cards and/or tubes)
- Up to 36 ID-Card and 24 tube capacity
- Panel PC touch-screen with easy-to-use interface (icon based)
- Multi-language software
- Batch specific timers (3 for ID-Cards, 2 for tubes)
- Standard program: 15 minutes at 37°C
- Programmable incubation time
- CE compliant according to IVD Directive 98/79/EC


ID-Card rack



Head 24 ID-Cards



Saxo centrifuge head



Ordering Information

Catalog #	Description
ID-Incubator L	
009203	ID-Incubator L, 1 unit